



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



MEDICAL EVACUATION (MEDEVAC)

...from point of injury to homeland hospital



4-FOLDED FIELD CATASTROPE STRETCHER

...a brief information

NODIN = Norwegian Development and Innovation



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



NATO MEDEVAC Stretcher (NT-820MVU)

NATO Stock No: NSN-6530-25-160-2914

In an Aeromedical Evacuation development project and design program, we search worldwide for a Field Stretcher meeting all requirements for aerial transport of patient, both civilian and military.

Since we were not able to find a stretcher meeting all requirements, we decided to develop a new field stretcher to be used in all kind of carriers; on Land, Sea and Air.

- Designed and tested iaw NATO Stanag 2040 and 3204
- Integrated 7 point quick release Harness
- 4 folded with telescopic handles
- Pinch free hinges
- Low storage volume
- 6 anchoring points for IV pole, arm-rest, medical bridge, light, hoist-straps etc.
- Minimum 10g crash load in all directions





NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



MEDEVAC Field Stretcher Mattress (NT-822MAD)

NATO Stock No: NSN-6530-25-161-8870



...unique mattress to reduce hypothermia risk and prevent bed sore for patient with limited circulation



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



MEDEVAC Field Stretcher Backrest (NT-850ABR)

NATO Stock No: NSN-6530-25-160-2915



Adjustable Backrest positions: 0°, 30°, 35°, 40° or 45°



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



MEDEVAC Field Stretcher and Mattress...

...stretcher alternatively packed in Bag or Backpack



...mattress packed in Dry Bag



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



NODIN IV-Pole (NT-821IVP) for NODIN MEDEVAC Stretcher



- 6 alternative and dedicated positions for mounting the IV-pole on the NODIN MEDEVAC Field Stretcher; 3 alternatives along each side of the stretcher
- Step less height adjustment from 410 mm to 740 mm
- Twist the IV-pole telescope left to release height position and twist right to lock in selected position.



NODIN AVIATION AS

YOUR NEEDS - OUR KNOWLEDGE - EQUALS INNOVATION



Nodin Aviation AS

PO Box 573, Nanset
N-3252 Larvik, NORWAY
NATO - NCAGE N4597

Phone: + 47 93048779

Email: post@medevac.no

www.medevac.no