

TECHNICAL DATA SHEET



Product: NODIN MEDEVAC STRETCHER (NT-820MVU)

NSN-6530-25-160-2914

Characteristics

Measurements according to NATO STANAG 2040
Load Capacity according to NATO STANAG 3204 AMD
Designed to meet Aerial Evacuation requirements
Pinch free hinges

Measurements

Total length , handles extended	2290mm
Total length , handles retracted	1980
Maximum width	584
Ground clearance	40

Folded stretcher

Height	520mm
Width	270mm
Depth	220mm

Weight:

9,7kg



Fabric:

Burning behavior:
Color

Self extinguishing. Passes FAR/CS 25.853a and ISO 3795:1999
Black

Patient Securing Harness


Burning behavior:
Color:
Buckles:

Self extinguishing. Passes FAR/CS 25.853a and ISO 3795:1999
Black
Female buckles: Aluminum/steel, male buckle stainless steel.

Stretcher frame:

Material
Color

Aluminum
Powder coated green



NODIN Aviation AS
Postbox 573, Nanset
3252 Larvik

www.medevac.no
post@medevac.no
phone: +47 41327943