

NT 822MAD

The NODIN Field stretcher mattress is a unique mattress designed to reduce the risk for hypothermia and prevent bed sores for patients with limited circulation.

- Integrated pocket for pillow
- Four carrying straps on each side
- Holes for harness
- Integrated straps for packing
- Straps for securing to stretcher handle

Secure the mattress to the stretcher by web-strap loop on the telescopic stretcher handle

Use an inflatable disposable pillow or other soft object to give the patient a more comfortable transport

The mattress is packed in a waterproof bag.

Height when packed:

~ 60 cm / 23.6"

Diameter when packed:

~ 25 cm / 9.8"

Weight of bag:

0.5 kg / 1.1 lb



Field Stretcher Mattress (NSN: 6530-25-161-8870)



Length: 185 cm - **Width:** 54 cm - **Thickness:** 3 cm - **Weight:** 2,5 kg.
Designed for the NODIN medevac Stretcher (NT-820MVU).
Compatible with stretchers designed iaw STANAG 2040



NODIN Aviation AS
Vangsveien 10
1814 ASKIM
NORWAY

Phone: +47 33 32 79 43
Email: post@medevac.no
Web: www.medevac.no
NCAGE: N4597



Org.no 879 913 802