

NT 827SEB

Stretcher Equipment Bridge

The NODIN Stretcher Equipment bridge is an add-on for the NODIN Medevac Stretcher enabling the user to mount medical equipment to the stretcher for treatment or monitoring of the patient.

It consists of a table, 3 poles and 5 straps

Ease of use

The equipment can be strapped to the equipment table with the included straps or attached to the poles with clamps.

5 straps (NT-828ESS) are included with the equipment bridge.

Dimensions when stored:


Length: 50 cm

Width: 50 cm

Depth: 5 cm

Weight: 3,5 kgs

Made from powder coated aluminium.



The Stretcher Equipment Bridge is a welcome add-on to any field examination installation or to secure equipment for transport along with the patient.

The poles on the Equipment Bridge can be further extended with the NT-837MEX and used for equipment rails like the NT-831MRS shown in the picture above, or attachment for IV-bags NT-834IVH.



NODIN Aviation AS
Vangsvæien 10
1814 ASKIM
NORWAY

Phone: +47 33 32 79 43
Email: post@medevac.no
Web: www.medevac.no
NCAGE: N4597



Org.no 879 913 802