



NT 990SWC

Stretcher Wheel Cart



The wheeled litter carrier is a lightweight casualty transport system that allows one person to transport a casualty strapped to a stretcher, and also performs as a litter stand in the battlefield or medical stations.

The revolutionary design is foldable, and can be totally disassembled for compact storage, easy shipping and transportation. The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern combat litters and is equipped with non-puncture wheels. A single user can operate the Stretcher Wheel Cart/Wheeled Litter Carrier. The Stretcher Carrier is equipped with clamps, which are easy to operate. The clamps secure and release the stretcher frame in seconds.



Technical Information

| | |
|---------------------------|--|
| Dimension (Deploy): | L143 x W70 x H92 (cm) |
| Dimension (Half Folded): | L126 x W70 x H50 (cm) |
| Dimension (Fully Folded): | L107 x W70 x H25 (cm) |
| Wheel dimensions | 50cm diameter |
| Weight: | 17.5kgs |
| Material type: | Manufactured with Altrunium® lightweight powder coated aluminium |

Manufactured by:

Nodin Aviation partner:
Golden Season Pte Ltd



NODIN Aviation AS
Vangsveien 10
1814 ASKIM
NORWAY

Phone: +47 33 32 79 43
Email: post@medevac.no
Web: www.medevac.no
NCAGE: N4597



Org.no 879 913 802